

## ATP4B antibody

### Product Information

Catalog No.:	FNab00701
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Required for stabilization and maturation of the catalytic proton pump alpha subunit and may also involved in cell adhesion and establishing epithelial cell polarity.

### Immunogen information

Immunogen:	ATPase, H <sup>+</sup> /K <sup>+</sup> exchanging, beta polypeptide
Synonyms:	Potassium-transporting ATPase subunit beta Gastric H(+)/K(+) ATPase subunit beta Proton pump beta chain ATP4B
Observed MW:	35 kDa
Uniprot ID :	P51164

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



SW620 cells were subjected to SDS PAGE followed by western blot with FNab00701(ATP4B antibody) at dilution of 1:1000