

# ATPAF1 antibody

#### **Product Information**

Catalog No.: FNab00730

Size: 100µg
Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months(Avoid repeated freeze / thaw cycles.)

# **Background**

May play an essential role for the assembly of the mitochondrial F1-F0 complex.

## Immunogen information

Immunogen: ATP synthase mitochondrial F1 complex assembly factor 1

Synonyms: ATP synthase mitochondrial F1 complex assembly factor 1|ATP11

homolog|ATPAF1|ATP11

Observed MW: 36 kDa
Uniprot ID: Q5TC12

## **Application**


Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:





mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00730(ATPAF1 antibody) at dilution of 1:500