

AXIN2 antibody

Product Information

Catalog No.:	FNab00753
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B(By similarity).

Immunogen information

Immunogen:	axin 2
Synonyms:	Axin-2 Axin-like protein (Axil) Axis inhibition protein 2 Conductin AXIN2
Observed MW:	93-100 kDa
Uniprot ID :	Q9Y2T1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Rat brain were subjected to SDS PAGE followed by western blot with FNab00753(AXIN2 Antibody) at dilution of 1:1000