

NPM1 antibody

Product Information

Catalog No.:	FNab00759
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nucleophosmin(nucleolar phosphoprotein B23, numatrin)
Synonyms:	Nucleophosmin (NPM) Nucleolar phosphoprotein B23 Nucleolar protein NO38 Numatrin NPM1 NPM
Observed MW:	38 kDa
Uniprot ID :	P06748

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:500-1:3000; IF: 1:50-1:500; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab00759(NPM1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00759(B23 Antibody) at dilution of 1:1000