

## B3GLCT antibody

### Product Information

Catalog No.:	FNab00767
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	beta 1,3-galactosyltransferase-like
Synonyms:	Beta-1,3-glucosyltransferase (Beta3Glc-T) Beta 3-glucosyltransferase Beta-3-glycosyltransferase-like B3GLCT B3GALTL B3GTL
Observed MW:	57 kDa
Uniprot ID :	Q6Y288

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human endometrial cancer using FNab00767(B3GALTL antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab00767(B3GALTL antibody) at dilution of 1:600