

BAALC antibody

Product Information

Catalog No.:	FNab00779
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a synaptic role at the postsynaptic lipid rafts by interacting with CAMK2A.

Immunogen information


Immunogen:	brain and acute leukemia, cytoplasmic
Synonyms:	Brain and acute leukemia cytoplasmic protein BAALC
Observed MW:	19-25 kDa
Uniprot ID :	Q8WXS3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00779(BAALC Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00779(BAALC Antibody) at dilution of 1:600