

ADGRB2-Specific antibody

Product Information

Catalog No.:	FNab00791
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


BAI2 belongs to the G-protein coupled receptor 2 family and LN-TM7 subfamily. BAI2 might be involved in angiogenesis inhibition. This antibody is specific to BAI2.

Immunogen information

Immunogen:	brain-specific angiogenesis inhibitor 2
Synonyms:	Adhesion G protein-coupled receptor B2 Brain-specific angiogenesis inhibitor 2 ADGRB2 BAI2
Observed MW:	None
Uniprot ID :	O60241

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab00791(BAI2-Specific Antibody) at dilution of 1:50