

## BAIAP2L1 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab00794  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


May function as adapter protein. Involved in the formation of clusters of actin bundles. Plays a role in the reorganization of the actin cytoskeleton in response to bacterial infection.

### Immunogen information


|              |  |
|--------------|--|
| Immunogen:   | BAI1-associated protein 2-like 1   |
| Synonyms:    | Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 1 (BAI1-associated protein 2-like protein 1) Insulin receptor tyrosine kinase substrate BAIAP2L1 IRTKS |
| Observed MW: | 57 kDa   |
| Uniprot ID : | Q9UHR4   |

### Application

|                       |                                       |
|-----------------------|---------------------------------------|
| Reactivity:           | Human                                 |
| Tested Application:   | ELISA, WB, IHC                        |
| Recommended dilution: | WB: 1:500 - 1:2000; IHC: 1:50 - 1:200 |
| Image:                |                                       |



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab00794( BAIAP2L1 antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00794( BAIAP2L1 Antibody) at dilution of 1:600