

BCAP29 antibody

Product Information

Catalog No.:	FNab00801
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	B-cell receptor-associated protein 29
Synonyms:	B-cell receptor-associated protein 29 (BCR-associated protein 29, Bap29) BCAP29 BAP29
Observed MW:	28 kDa
Uniprot ID :	Q9UHQ4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab00801(BCAP29 antibody) at dilution of 1:200



Raji cells were subjected to SDS PAGE followed by western blot with FNab00801(BCAP29 antibody) at dilution of 1:1000