

BARHL2 antibody

Product Information

Catalog No.:	FNab00804
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potential regulator of neural basic helix-loop-helix genes.

Immunogen information


Immunogen:	BarH-like homeobox 2
Synonyms:	BarH-like 2 homeobox protein BARHL2
Observed MW:	42-45 kDa
Uniprot ID :	Q9NY43

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00804 (BARHL2 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00804 (BARHL2 Antibody) at dilution of 1:300