

BCL2L10 antibody

Product Information

Catalog No.:	FNab00844
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Promotes cell survival. Suppresses apoptosis induced by BAX but not BAK.

Immunogen information


Immunogen:	BCL2-like 10(apoptosis facilitator)
Synonyms:	Bcl-2-like protein 10 (Bcl2-L-10) Anti-apoptotic protein Boo Anti-apoptotic protein NrH Apoptosis regulator Bcl-B BCL2L10 BCL-B BCLB BOO DIVA
Observed MW:	22 kDa
Uniprot ID :	Q9HD36

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung using FNaB00844(BCL2L10 antibody) at dilution of 1:50



A2780 cells were subjected to SDS PAGE followed by western blot with FNaB00844(BCL2L10 antibody) at dilution of 1:200