

BCL2L12 antibody

Product Information

Catalog No.:	FNab00845
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	BCL2-like 12(proline rich)
Synonyms:	Bcl-2-like protein 12 (Bcl2-L-12) Bcl-2-related proline-rich protein BCL2L12 BPR
Observed MW:	32 kDa
Uniprot ID :	Q9HB09

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer slide using FNab00845(BCL2L12 Antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab00845(BCL2L12 Antibody) at dilution of 1:1000