

BCO1 antibody

Product Information

Catalog No.:	FNab00851
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	beta-carotene 15,15'-monooxygenase 1
Synonyms:	Beta,beta-carotene 15,15'-dioxygenase Beta-carotene dioxygenase 1 Beta-carotene oxygenase 1 BCO1 BCDO BCDO1 BCMO1
Observed MW:	65 kDa
Uniprot ID :	Q9HAY6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00851(BCMO1 Antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab00851(BCMO1 Antibody) at dilution of 1:300