

ABCG2 antibody

Product Information

Catalog No.:	FNab00855
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ATP-binding cassette, sub-family G(WHITE), member 2
Synonyms:	Broad substrate specificity ATP-binding cassette transporter ABCG2 ATP-binding cassette sub-family G member 2 Breast cancer resistance protein CDw338 Mitoxantrone resistance-associated protein Placenta-specific ATP-binding cassette transporter Urate exporter ABCG2 ABCP BCRP BCRP1 MXR
Observed MW:	60-70 kDa
Uniprot ID :	Q9UNQ0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab00855(BCRP,ABCG2 Antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00855(BCRP,ABCG2 antibody) at dilution of 1:300