

BEND4 antibody

Product Information

Catalog No.:	FNab00862
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	BEN domain containing 4
Synonyms:	BEN domain-containing protein 4 Coiled-coil domain-containing protein 4 BEND4 CCDC4
Observed MW:	55-58 kDa;33-36 kDa
Uniprot ID :	Q6ZU67

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab00862(BEND4 Antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00862(BEND4 Antibody) at dilution of 1:1000