

BEST3 antibody

Product Information

Catalog No.:	FNab00867
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	bestrophin 3
Synonyms:	Bestrophin-3 Vitelliform macular dystrophy 2-like protein 3 BEST3 VMD2L3
Observed MW:	76-78 kDa, 50-55 kDa
Uniprot ID :	Q8N1M1

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using FNab00867(Bestrophin-3 antibody at dilution of 1:200



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab00867(Bestrophin-3 antibody) at dilution of 1:800