

## SGCB antibody

### Product Information

Catalog No.:	FNab00878
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	sarcoglycan, beta(43kDa dystrophin-associated glycoprotein)
Synonyms:	Beta-sarcoglycan (Beta-SG) 43 kDa dystrophin-associated glycoprotein (43DAG) A3b SGCB
Observed MW:	48 kDa
Uniprot ID :	Q16585

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00878(Beta sarcoglycan antibody) at dilution of 1:600