

## BHMT2 antibody

### Product Information

Catalog No.:	FNab00891
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the regulation of homocysteine metabolism. Converts homocysteine to methionine using S-methylmethionine(SMM) as a methyl donor.

### Immunogen information

Immunogen:	betaine-homocysteine methyltransferase 2
Synonyms:	S-methylmethionine--homocysteine S-methyltransferase BHMT2 (SMM-hcy methyltransferase) Betaine--homocysteine S-methyltransferase 2 BHMT2
Observed MW:	40-45 kDa
Uniprot ID :	Q9H2M3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:3000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00890(BHMT2 antibody) at dilution of 1:2000