

## BID antibody

### Product Information

Catalog No.:	FNab00893
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6E8
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


The major proteolytic product p15 BID allows the release of cytochrome c(By similarity). Isoform 1, isoform 2 and isoform 4 induce ICE-like proteases and apoptosis. Isoform 3 does not induce apoptosis. Counters the protective effect of Bcl-2.

### Immunogen information

Immunogen:	BH3 interacting domain death agonist
Synonyms:	BH3-interacting domain death agonist p22 BID (BID) BH3-interacting domain death agonist p15 Alternative names: p15 BID BH3-interacting domain death agonist p13 Alternative names: p13 BID BH3-interacting domain death agonist p11 Alternative names: p11 BID BID
Observed MW:	22 kDa
Uniprot ID :	P55957

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:2000-1:10000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00893( BID Antibody) at dilution of 1:3000