

CEL antibody

Product Information

Catalog No.:	FNab00896
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	carboxyl ester lipase(bile salt-stimulated lipase)
Synonyms:	Bile salt-activated lipase (BAL) Bile salt-stimulated lipase (BSSL) Bucelipase Carboxyl ester lipase Cholesterol esterase Pancreatic lysophospholipase Sterol esterase CEL BAL
Observed MW:	67 kDa
Uniprot ID :	P19835

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human liver tissue were subjected to SDS PAGE followed by western blot with FNab00896(Bile salt-activated lipase antibody) at dilution of 1:500