

## BIRC8 antibody

### Product Information

Catalog No.:	FNab00901
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Protects against apoptosis mediated by BAX.

### Immunogen information

Immunogen:	baculoviral IAP repeat-containing 8
Synonyms:	Baculoviral IAP repeat-containing protein 8 Inhibitor of apoptosis-like protein 2 (IAP-like protein 2, ILP-2) Testis-specific inhibitor of apoptosis BIRC8 ILP2
Observed MW:	27 kDa
Uniprot ID :	Q96P09

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00901(BIRC8 antibody) at dilution of 1:1000