

## GDF2 antibody

### Product Information

Catalog No.:	FNab00921
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	growth differentiation factor 2
Synonyms:	Growth/differentiation factor 2 (GDF-2) Bone morphogenetic protein 9 (BMP-9) GDF2 BMP9
Observed MW:	30 kDa, 46 kDa
Uniprot ID :	Q9UK05

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using FNab00921(BMP9 Antibody) at dilution of 1:50



L02 cells were subjected to SDS PAGE followed by western blot with FNab00921(BMP9 Antibody) at dilution of 1:300