

UBL7 antibody

Product Information

Catalog No.:	FNab00925
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ubiquitin-like 7(bone marrow stromal cell-derived)
Synonyms:	Ubiquitin-like protein 7 Bone marrow stromal cell ubiquitin-like protein (BMSC-UbP) Ubiquitin-like protein SB132 UBL7 BMSCUBP
Observed MW:	41 kDa
Uniprot ID :	Q96S82

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab00925(UBL7 antibody) at dilution of 1:50



IP Result of anti-BMSC UbP (IP:FNab00925, 4ug; Detection:FNab00925 1:1000) with HeLa cells lysate 2000ug.



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00925(UBL7 antibody) at dilution of 1:400