

## TBXT antibody

### Product Information

Catalog No.:	FNab00940
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	T, brachyury homolog(Mouse)
Synonyms:	T-box transcription factor T Brachyury protein Protein T TBXT T
Observed MW:	47 kDa
Uniprot ID :	O15178

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	



H460 cells were subjected to SDS PAGE followed by western blot with FNab00940(Brachyury-T antibody) at dilution of 1:500