

## BRF2 antibody

### Product Information

Catalog No.:	FNab00955
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


General activator of RNA polymerase III transcription. Factor exclusively required for RNA polymerase III transcription of genes with promoter elements upstream of the initiation sites. It recognizes a 50-60kd band in WB analysis.

### Immunogen information

Immunogen:	BRF2, subunit of RNA polymerase III transcription initiation factor, BRF1-like
Synonyms:	Transcription factor IIIB 50 kDa subunit (TFIIIB50, hTFIIIB50) B-related factor 2 (BRF-2) hBRFU BRF2 BRFU
Observed MW:	50-60 kDa
Uniprot ID :	Q9HAW0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



HT-1080 cells were subjected to SDS PAGE followed by western blot with FNab00955(BRF2 antibody) at dilution of 1:500