

PTK6 antibody

Product Information

Catalog No.:	FNab00959
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	PTK6 protein tyrosine kinase 6
Synonyms:	Protein-tyrosine kinase 6 Breast tumor kinase Tyrosine-protein kinase BRK PTK6 BRK
Observed MW:	52 kDa
Uniprot ID :	Q13882

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab00959(PTK6 antibody) at dilution of 1:100



IP Result of anti-BRK (IP:FNab00959, 4ug; Detection:FNab00959 1:800) with HeLa cells lysate 3440ug.



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00959(PTK6 antibody) at dilution of 1:500