

BRMS1L antibody

Product Information

Catalog No.:	FNab00961
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the histone deacetylase(HDAC1)-dependent transcriptional repression activity. When overexpressed in lung cancer cell line that lacks p53/TP53 expression, inhibits cell growth.

Immunogen information

Immunogen:	breast cancer metastasis-suppressor 1-like
Synonyms:	Breast cancer metastasis-suppressor 1-like protein BRMS1-homolog protein p40 BRMS1-like protein p40 BRMS1L
Observed MW:	40 kDa
Uniprot ID :	Q5PSV4

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab00961(BRMS1L Antibody) at dilution of 1:600