

POU3F2-Specific antibody

Product Information

Catalog No.:	FNab00963
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	POU class 3 homeobox 2
Synonyms:	POU domain, class 3, transcription factor 2 Brain-specific homeobox/POU domain protein 2 (Brain-2, Brn-2) Nervous system-specific octamer-binding transcription factor N-Oct-3 Octamer-binding protein 7 (Oct-7) Octamer-binding transcription factor 7 (OTF-7) POU3F2 BRN2 OCT7 OTF7
Observed MW:	47 kDa
Uniprot ID :	P20265

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab00963(POU3F2-Specific antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00963(POU3F2-Specific antibody) at dilution of 1:600