

HGH1 antibody

Product Information

Catalog No.:	FNab00964
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 8 open reading frame 30A
Synonyms:	Protein HGH1 homolog HGH1 C8orf30A C8orf30B FAM203A FAM203B BRP16
Observed MW:	42-45 kDa
Uniprot ID :	Q9BTY7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with FNab00964(BRP16 Antibody) at dilution of 1:300