

BTBD1 antibody

Product Information

Catalog No.:	FNab00971
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable substrate-specific adapter of an E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Immunogen information


Immunogen:	BTB(POZ) domain containing 1
Synonyms:	BTB/POZ domain-containing protein 1 Hepatitis C virus NS5A-transactivated protein 8 (HCV NS5A-transactivated protein 8) BTBD1 C15orf1 NS5ATP8
Observed MW:	53 kDa
Uniprot ID :	Q9H0C5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab00971(BTBD1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00971(BTBD1 antibody) at dilution of 1:500