

BTN2A1 antibody

Product Information

Catalog No.:	FNab00981
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	butyrophilin, subfamily 2, member A1
Synonyms:	Butyrophilin subfamily 2 member A1 BTN2A1 BT2.1 BTF1
Observed MW:	69 kDa
Uniprot ID :	Q7KYR7

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00981(BTN2A1 antibody) at dilution of 1:600