

BTN2A2 antibody

Product Information

Catalog No.:	FNab00982
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibits the proliferation of CD4 and CD8 T-cells activated by anti-CD3 antibodies, T-cell metabolism and IL2 and IFNG secretion.

Immunogen information

Immunogen:	butyrophilin, subfamily 2, member A2
Synonyms:	Butyrophilin subfamily 2 member A2 BTN2A2 BT2.2 BTF2
Observed MW:	59 kDa
Uniprot ID :	Q8WVV5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00982(BTN2A2 antibody) at dilution of 1:500