

BTNL2 antibody

Product Information

Catalog No.:	FNab00985
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Negative regulator of T-cell proliferation.

Immunogen information

Immunogen:	butyrophilin-like 2
Synonyms:	Butyrophilin-like protein 2 (BTL-II) BTNL2
Observed MW:	65 kDa
Uniprot ID :	Q9UIR0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00985(BTNL2 Antibody) at dilution of 1:600