

BTNL8 antibody

Product Information

Catalog No.:	FNab00986
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May stimulate primary immune response. Acts on T-cell stimulated sub-optimally through the TCR/CD3 complex stimulating their proliferation and cytokine production.

Immunogen information

Immunogen:	butyrophilin-like 8
Synonyms:	Butyrophilin-like protein 8 BTNL8
Observed MW:	50 kDa
Uniprot ID :	Q6UX41

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab00986(BTNL8 Antibody) at dilution of 1:600