

BZW1 antibody

Product Information

Catalog No.:	FNab00995
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Enhances histone H4 gene transcription but does not seem to bind DNA directly.

Immunogen information


Immunogen:	basic leucine zipper and W2 domains 1
Synonyms:	eIF5-mimic protein 2 Basic leucine zipper and W2 domain-containing protein 1 Protein OrfBZW1 5MP2 BZAP45 KIAA0005
Observed MW:	48 kDa
Uniprot ID :	Q7L1Q6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab00995(BZW1 antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00995(BZW1 antibody) at dilution of 1:500