

## TBATA antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab00997  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information


|              |  |
|--------------|--|
| Immunogen:   | chromosome 10 open reading frame 27  |
| Synonyms:    | Protein TBATA Protein SPATIAL Stromal protein associated with thymii and lymph node homolog Thymus, brain and testes-associated protein TBATA C10orf27 SPATIAL |
| Observed MW: | 39 kDa   |
| Uniprot ID : | Q96M53   |

### Application

|                       |                                   |
|-----------------------|-----------------------------------|
| Reactivity:           | Human                             |
| Tested Application:   | ELISA, WB, IHC                    |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200 |
| Image:                |                                   |



Immunohistochemistry of paraffin-embedded human testis slide using FNab00997(C10orf27 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00997(C10orf27 Antibody) at dilution of 1:600