

## C10orf67 antibody

### Product Information

Catalog No.:	FNab01000
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	chromosome 10 open reading frame 67
Synonyms:	Uncharacterized protein C10orf67, mitochondrial C10orf67 LINC01552
Observed MW:	20 kDa
Uniprot ID :	Q8IYJ2

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:1000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01000( C10orf67 Antibody) at dilution of 1:300