

ANAPC15 antibody

Product Information

Catalog No.:	FNab01004
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 11 open reading frame 51
Synonyms:	Anaphase-promoting complex subunit 15 (APC15) ANAPC15 C11orf51
Observed MW:	22 kDa
Uniprot ID :	P60006

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung using FNab01004(C11orf51 antibody) at dilution of 1:50



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab01004(C11orf51 antibody) at dilution of 1:300