

## AAMDC antibody

### Product Information

Catalog No.:	FNab01007
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 11 open reading frame 67
Synonyms:	Mth938 domain-containing protein Adipogenesis associated Mth938 domain-containing protein AAMDC C11orf67
Observed MW:	13 kDa
Uniprot ID :	Q9H7C9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary using FNab01007(C11orf67 Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab01007 (C11orf67 antibody) at dilution of 1:1000