

## IFTAP antibody

### Product Information

Catalog No.:	FNab01009
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None


### Immunogen information

Immunogen:	chromosome 11 open reading frame 74
Synonyms:	Intraflagellar transport-associated protein Protein HEPIS IFTAP C11orf74 NWC
Observed MW:	25 kDa
Uniprot ID :	Q86VG3


### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	





Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab01009(C11orf74 Antibody) at dilution of 1:50



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab01009(C11orf74 Antibody) at dilution of 1:600