

## SDHAF2 antibody

### Product Information

Catalog No.:	FNab01010
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 11 open reading frame 79
Synonyms:	Succinate dehydrogenase assembly factor 2, mitochondrial (SDH assembly factor 2, SDHAF2) SDHAF2 C11orf79 PGL2 SDH5
Observed MW:	30 kDa
Uniprot ID :	Q9NX18

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab01010(C11orf79 Antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab01010(C11orf79 Antibody) at dilution of 1:300