

## C12orf54 antibody

### Product Information

Catalog No.:	FNab01013
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	chromosome 12 open reading frame 54
Synonyms:	Uncharacterized protein C12orf54 C12orf54
Observed MW:	14 kDa
Uniprot ID :	Q6X4T0

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue using FNab01013(C12orf54 antibody) at dilution of 1:50