

## KATNBL1 antibody

### Product Information

Catalog No.:	FNab01028
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 15 open reading frame 29
Synonyms:	KATNB1-like protein 1 Katanin p80 subunit B-like 1 KATNBL1 C15orf29
Observed MW:	28 kDa
Uniprot ID :	Q9H079

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab01028(C15orf29 Antibody) at dilution of 1:50



Neuro-2a cells were subjected to SDS PAGE followed by western blot with FNab01028(C15orf29 Antibody) at dilution of 1:600