

CMC2 antibody

Product Information

Catalog No.:	FNab01035
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None


Immunogen information

Immunogen:	chromosome 16 open reading frame 61
Synonyms:	COX assembly mitochondrial protein 2 homolog CMC2 C16orf61
Observed MW:	9 kDa
Uniprot ID :	Q9NRP2

Application

Reactivity:HumanTested Application:ELISA, WBRecommended dilution:WB: 1:500-1:2000Image:Image:





PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01035(C16orf61 Antibody) at dilution of 1:600