

HAPSTR1 antibody

Product Information

Catalog No.:	FNab01037
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 16 open reading frame 72
Synonyms:	HUWE1-associated protein modifying stress responses 1 Telomere attrition and p53 response 1 protein HAPSTR1 C16orf72 TAPR1
Observed MW:	31 kDa
Uniprot ID :	Q14CZ0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01037(C16orf72 Antibody) at dilution of 1:300