

## CTU2 antibody

### Product Information

Catalog No.:	FNab01038
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 16 open reading frame 84
Synonyms:	Cytoplasmic tRNA 2-thiolation protein 2 Cytosolic thiouridylase subunit 2 CTU2 C16orf84 NCS2
Observed MW:	66 kDa
Uniprot ID :	Q2VPK5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab01038(C16orf84 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab01038(C16orf84 Antibody) at dilution of 1:600