

ELP5 antibody

Product Information

Catalog No.:	FNab01041
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	chromosome 17 open reading frame 81
Synonyms:	Elongator complex protein 5 Dermal papilla-derived protein 6 S-phase 2 protein ELP5 C17orf81 DERP6
Observed MW:	None
Uniprot ID :	Q8TE02

Application

Reactivity:	Human
Tested Application:	ELISA
Recommended dilution:	None