

TIMM21 antibody

Product Information

Catalog No.:	FNab01045
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 18 open reading frame 55
Synonyms:	Mitochondrial import inner membrane translocase subunit Tim21 TIM21-like protein, mitochondrial TIMM21 C18orf55 TIM21
Observed MW:	None
Uniprot ID :	Q9BVV7

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human liver using FNab01045(C18orf55 antibody) at dilution of 1:100