

COA7 antibody

Product Information

Catalog No.:	FNab01062
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 1 open reading frame 163
Synonyms:	Cytochrome c oxidase assembly factor 7 Beta-lactamase hcp-like protein Respiratory chain assembly factor 1 Sel1 repeat-containing protein 1 COA7 C1orf163 RESA1 SELRC1
Observed MW:	25-28 kDa
Uniprot ID :	Q96BR5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01062(C1orf163 Antibody) at dilution of 1:1000