

ODR4 antibody

Product Information

Catalog No.:	FNab01064
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 1 open reading frame 27
Synonyms:	Protein odr-4 homolog (hODR-4) LAG1-interacting protein Transactivated by transforming growth factor beta protein 1 ODR4 C1orf27 TTG1 TTG1A
Observed MW:	50 kDa
Uniprot ID :	Q5SWX8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01064(C1orf27 Antibody) at dilution of 1:300